



Test kit, ozone, model OZ-2

Kat. číslo produktu:

2064400

CZK Cena bez DPH:

Kontaktujte nás

Dodání do 2 týdnů

Nebezpečné



Hach's most economical test kit for ozone measurement utilizes the simple, low-cost color disc method.

Hach's most economical test kit for ozone measurement

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Technické údaje

Metoda:	Color disc/DPD
Parameter:	Ozone
Parametr:	Ozone - As O ₃
Platforma:	Barevný disk
Počet testů:	100
Přístroj:	OZ-2
Rozsah měření:	0 - 2,2 mg/L O ₃
Smallest increment:	0.05
Standardní rozsah měření:	0 - 2.2 mg/L O ₃
Styl kufříku:	B
Typ testovací soupravy:	OZ-2