



## Beaker, polypropylene, low form, 400 mL, 6/pk

|                          |                 |
|--------------------------|-----------------|
| Kat. číslo produktu:     | 108078          |
| CZK Cena bez DPH:        | Kontaktujte nás |
| Datum odeslání neuvedeno |                 |

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes  $\geq 250$  mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

---

### Technické údaje

|                  |                         |
|------------------|-------------------------|
| Jednotky:        | 6 /pk                   |
| Kapacita:        | 400 mL                  |
| Materiál:        | Polypropylene, Low Form |
| Materiál skříně: | Polypropylene           |
| Průměr:          | 99.1 mm with spout      |