



## Stablcal Stabilizovaný primární formazínový standard pro měření turbidity, 20 NTU

**Kat. číslo produktu:**

**2660149**

**Nebezpečné**

CZK Cena bez DPH:  
Dostupné

Kontaktujte nás



Připraveno k použití, stabilizovaný primární formazínový standard v 500 ml lahvi, 20 NTU.

### Precise, Pre-mixed Formazin Solutions

With proprietary manufacturing technology, Hach prepares Stablcal Standards in precise concentrations at a wide variety of concentrations - providing accurate turbidimeter calibrations.

### Based on Formazin Primary Standards

Stablcal Stabilised Formazin Standards contain the same light scattering polymer as traditional Formazin primary turbidity standards - and perform equivalently.

### Enhanced Stability Saves Time and Labor

While Formazin primary standards can be prepared directly in the lab, the process is labor intensive and time consuming - and the diluted standards are not stable. Stablcal is a stabilised version of Formazin - providing repeatable results, without the time and effort required for lab-prepared standards.

### Confidence in results

Hach was the first to apply Formazin as a turbidity standard in the drinking water industry, more than 40 years ago. Stablcal standards are based on decades of experience in turbidity measurement.

---

## Technické údaje

Koncentrace: 20 NTU

Platforma: Roztok