



Amino Acid F Reagent for Silica analysis, 2.9 L

Kat. číslo produktu:

2353103

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Dostupné



Amino Acid F Reagent (2.7 L) used in the Series 5000 Silica Analyser.

Amino Acid F reagent is supplied in two parts to ensure a long shelf life: a 3-liter bottle of Amino Acid F Solution and a bottle of Amino Acid F Powder. Powder must be added to the solution immediately before installation in the analyser.

Technické údaje

Obsah balení:

1 x 2.7 L & 1 x 308 g

Parametr:

Silica

Podmínky skladování:

Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Přístroj:

5500sc Silica

Rozsah měření:

N/A