



Pall Membrane filter, pore size 0.45 µm, diameter 47 mm, 100/pk

Kat. číslo produktu:

1353000

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno



Ensure complete retention and maximum recovery of bacteria. Plain, non-sterile.

These 47 mm diameter, 0.45 µm pore size GN-Metricel membrane filters meet or surpass all U.S. Environmental Protection Agency (USEPA) and American Society for Testing and Materials (ASTM) specifications to ensure complete retention and maximum

recovery of bacteria. Black membranes for yeast and mold are also available. Black membranes make yeast and mold colonies more distinct and easier to count. Filters are available in convenient, presterilized, individually sealed packs (P/N 13530-01

and 24903-00) or in resealable envelopes with 10 filters and 10 absorbent pads each (P/N 21269-00). Phosphate free.

Technické údaje

Filter Type:	Filter
Materiál:	Mixed-Cellulose Ester
Obsah balení:	100 pk
Průměr:	47 mm
Velikost póru:	0.45 µm