



FILTER, MEMB GRID 47MM PK/200 GELMAN

Kat. číslo produktu: 1353001
CZK Cena bez DPH: Kontaktujte nás
Datum odeslání neuvedeno

Individually packaged, gridded, sterile, pk/200

These 47 mm diameter, 0.45 µm pore size GN-Metricel membrane filters meet or surpass all U.S. Environmental Protection Agency (USEPA) and American Society for Testing and Materials (ASTM) specifications to ensure complete retention and maximum recovery of bacteria. 200/pk, without pad. Gelman/Pall.

Black membranes for yeast and mold are also available. Black membranes make yeast and mold colonies more distinct and easier to count. Filters are available in convenient, presterilized, individually sealed packs (P/N 1353001 and 2490300) or in resealable envelopes with 10 filters and 10 absorbent pads each (P/N 2126900). Phosphate free.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of

your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach offers a variety of media and apparatus to meet all of

your recreational testing needs.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking

water samples, Hach has prepared media and convenient apparatus to make

testing easy and fast.

Technické údaje

| | |
|----------------|-----------------------|
| Filter Type: | Filter with grid |
| Materiál: | Mixed-Cellulose Ester |
| Obsah balení: | 200 pk |
| Průměr: | 47 mm |
| Velikost póru: | 0.45 µm |