



## ISE25Cl-9 Chloride ion selective electrode (ISE) half-cell, screw cap

<b>Kat. číslo produktu:</b>	<b>E41M003</b>
CZK Cena bez DPH:	Kontaktujte nás
Dodání do 2 týdnů	

### The right probe for your application.

To meet your analysis needs, the Radiometer Analytical ion selective electrode range offers high quality sensors for a variety of ionic species in solution.

### Factory controlled to ensure high quality and premium performance.

Each Radiometer Analytical Ion Selective Electrode is carefully controlled in production. Each ISE is delivered with an individual control report including information on the response curve, checkpoint value, and slope.

---

## Technické údaje

Záruka:	6 měsíců
---------	----------