



pH glass electrode, refillable, pure water applications

Kat. číslo produktu:

LZY025

CZK Cena bez DPH:

Kontaktujte nás

Dostupné

pH electrode for use in ultra-pure water

Specially developed for use in ultra-pure water

Refillable glass electrodes for special requirements and high-precision measurements

Improved durability through triple ceramic diaphragm

Easy to install in stainless steel or polypropylene flow fittings

Technické údaje

Čas odezvy :	< 30 s
Délka:	120 mm
Elektrolyt:	KCl liquid
Hmotnost:	0.15 kg
Koncovka kabelu:	Plug head S7
Materiál (rukojeť elektrody):	Glass
Membrána:	Triple ceramic
Provozní teplotní rozsah:	-30 - 80 °C
Reference system:	Ag/AgCl
Rozsah měření:	0 - 12 pH
Teplotní senzor:	None
Záruka:	6 měsíců