



Beaker, polypropylene, low form, 600 mL

Kat. číslo produktu:

108052

CZK Cena bez DPH:

Kontaktujte nás

Dodání do 2 týdnů

Griffin low form. For general laboratory use. Non-toxic, noncontaminating, and biologically inert. Autoclavable. Combines no drip pouring, durability, and translucency. Approximate volume scales molded in for small sizes. Sizes 250 mL and larger feature silk-screened graduations and double spout for easier pouring. Nalgene.

Technické údaje

Číslo metody:	8006
Jednotky:	1 each
Kapacita:	600 mL
Materiál:	Polypropylene, Low Form
Přístroj:	All except DR/4000, PC II
Průměr:	110.7 mm with spout