



## Beaker, polypropylene, low form, 1000 mL

|                      |                 |
|----------------------|-----------------|
| Kat. číslo produktu: | 108053          |
| CZK Cena bez DPH:    | Kontaktujte nás |
| Dodání do 2 týdnů    |                 |

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes  $\geq 250$  mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

---

### Technické údaje

|                  |                         |
|------------------|-------------------------|
| Jednotky:        | 1 each                  |
| Kapacita:        | 1000 mL                 |
| Materiál:        | Polypropylene, Low Form |
| Materiál skříně: | Polypropylene           |
| Průměr:          | 129.8 mm with spout     |