



**Be Right™**



## **Funnel, analytical, 254 mL approx. Volume**

**Kat. číslo produktu:**

**108370**

**CZK Cena bez DPH:**

**Kontaktujte nás**

**Datum odeslání neuvedeno**

Standard long stem funnel meets the needs of a wide variety of applications. Mounded ribbing prevents "air lock" commonly occurring during funneling and filtrations. A steep, 58-degree cone slope makes for fast and easy flow of either liquids or small solids. Polypropylene construction is resistant to chemicals and makes this funnel fully autoclavable for sterilization.

---

### **Technické údaje**

Obsah balení:

1 each