



Flask, volumetric, polypropylene, 100 mL

Kat. číslo produktu:	1406042
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

Polypropylene flask and screw cap. Durable with excellent chemical resistance, with the same accuracy as glass. Permanent color graduation line with marking spot. Calibrated “to contain/deliver” at 20°C. The 200/203 capacity flask has dual marks useful in DO determination via Winkler Titration.

Technické údaje

Certifikace:	Not Applicable
Jednotky:	1 each
Kapacita:	100 mL
Materiál skříně:	Polypropylene