



Test kit, nitrogen, nitrate color disc, model NI-14

Kat. číslo produktu:

1416100

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Dodání do 2 týdnů



Utilizing the simple color disc method, this kit measures Nitrogen, Nitrate as NO₃ and contains 50 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Metoda:	Color disc/ Cadmium reduction
Parameter:	Nitrogen, Nitrate
Parametr:	Nitrogen, Nitrate - As NO ₃ ⁻ -N
Platforma:	Barevný disk
Počet testů:	50
Přístroj:	NI-14
Rozsah měření:	Dual Range
Rozsah měření (2):	0 - 10 mg/L NO ₃ ⁻ -N
Smallest increment:	0.02 0.2
Standardní rozsah měření:	0 - 1 mg/L NO ₃ ⁻ -N; 0 - 10 mg/L NO ₃ ⁻ -N
Styl kufříku:	B
Typ testovací soupravy:	NI-14**