



Be Right™



m-FC broth, dehydrated, 100 g

Kat. číslo produktu:

1462426

ZASTARALÁ POLOŽKA

Tato položka již není dostupná.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of

your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make

wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking

water samples, Hach has prepared media and convenient apparatus to make

testing easy and fast.

Technické údaje

Media Type:

m-FC Broth

Obsah balení:

100 g