



Iron and manganese color disc test kit, model IR-20

Kat. číslo produktu:

146300

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno

Nebezpečné



When customers complain of stained and dingy laundry, this kit lets you do quick, on-the-spot tests for iron and manganese to see where the problem lies.

Method / Chemistry: Color disc colorimetry. Test: Iron/ Manganese. Range (mg/L): 0-5/ 0-3. Smallest I

Simple, fast, color disc method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Aplikace:	Water Quality
Nejmenší dílky přírustku:	0.2
Parametr:	Iron and Manganese
Platforma:	Barevný disk
Přístroj:	IR-20
Rozsah měření:	0 - 3,0
Rozsah měření (2):	0 - 4,0 mg/L Fe