



Be Right™



Sběrač stojanů 24 pozic

Kat. číslo produktu:

17012

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno

24 position no wire test tube rack. 25mm OD holes. rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Outside dimensions are all the same at 246 x 104 x 64 mm high (9.7 in x 4.1 in x 2.5 in). (Lachat)