



Chromium color disc test kit, model CH-8

Kat. číslo produktu:

183400

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno



Widely used as a corrosion inhibitor in cooling towers and engine cooling systems, chromium is also used in surface finishing, leather tanning, and other industries. The most widely used form, hexavalent chromium, is very toxic and strict limits are placed on discharges. Utilizing the color disc method to determine chromium (Range: 0.1-1.5 mg/L), this kit contains 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Metoda:	Color disc/ Diphenylcarbazide
Parameter:	Chromium Hexavalent
Parametr:	Chromium, low range - As hexavalent Cr
Platforma:	Barevný disk
Počet testů:	100
Přístroj:	CH-8
Rozsah měření:	0 - 1,3 Cr ⁶⁺
Smallest increment:	0.1
Styl kufříku:	D
Typ testovací soupravy:	CH-8