



Cyanuric acid test kit, model CY-3

Kat. číslo produktu: 185102
CZK Cena bez DPH: Kontaktujte nás
Datum odeslání neuvedeno

A stabilizing agent for chlorine, cyanuric acid is used in swimming pool treatment. Cyanuric Acid Test Kit measures the turbidity formed in proportion to Cyanuric Acid concentration, after the reagent is added. Contains 100 tests.

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Aplikace:	Pool and Spa
Metoda:	Turbidimetry
Parametr:	Cyanuric Acid
Počet testů:	100
Přístroj:	CY-3
Rozsah měření:	20 - 100 mg/L
Smallest increment:	20
Styl kufříku:	D
Typ testovací soupravy:	CY-3