



Be Right™



Test kit tannin/lignin

Kat. číslo produktu:

193701

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno

Nebezpečné



A product of decomposed plant material, tannin occurs in natural waters and is used for boiler and cooling tower water treatment, dyeing, ink manufacture, tanning, and paper sizing. Lignin, also a plant product, is found in paper pulp manufacturing waste. Tannin and lignin are analyzed together and reported as mg/L tannic acid. This kit utilizes the simple, low-cost color disc method - the most accurate visual test kit method.

Tannin and lignin are analyzed together and reported as mg/L tannic acid.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Technické údaje

Metoda:	Color disc/Tyrosine
Parameter:	Tannin/Lignin
Parametr:	Tannin/Lignin - As tannic acid
Platforma:	Barevný disk
Počet testů:	100
Přístroj:	Tannin/Lignin
Rozsah měření:	0 - 15 mg/L as Tannic Acid
Rozsah měření (2):	0 - 150 mg/L
Smallest increment:	0.5
	5
Standardní rozsah měření:	0 - 15 mg/L as Tannic Acid; 0 - 150 mg/L
Styl kufříku:	R
Typ testovací soupravy:	TA-3