



Be Right™



Bottle, Sampling, High Density, Polyethylene, 2 L

Kat. číslo produktu:	2087154
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

High-density polyethylene bottle with wide-mouth (89 mm) opening. Excellent sampling bottle for grab samples, storing dry chemicals, and shipping liquids.

Technické údaje

Body Material:	High Density Plastic
Kapacita:	2 L
Materiál:	HDPE
Obsah balení:	1 each
Package size UOM:	1pk