



Be Right™



Flask, Erlenmeyer, glass w/screw cap, 125 mL

Kat. číslo produktu:

2089643

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

Technické údaje

Jednotky:

1 each

Kapacita:

125 mL

Materiál:

Glass Screw Cap