



Be Right™



Flask, Erlenmeyer, glass w/screw cap, 250 mL

Kat. číslo produktu:

2089646

ZASTARALÁ POLOŽKA

Tato položka již není dostupná.

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

Technické údaje

Jednotky:	1 each
Kapacita:	250 mL
Materiál:	Glass Screw Cap