



Flask, Erlenmeyer, polymethylpentene capacity 125 mL, 6/pk

Kat. číslo produktu:	2089873
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

Polymethylpentene flask with polypropylene screwcap features clarity, chemical resistance, strength, and autoclavability (Do not engage threads during cycle). Useful for titrations, preparation and storage of media.

Technické údaje

Jednotky:	6/pk
Kapacita:	125 mL
Materiál:	Polymethylpentene
Typ:	Erlenmeyer