



Be Right™



Centrifuge Tube, Conical Bottom, 10/pk

Kat. číslo produktu:	2115839
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

Centrifuge tube with conical bottom. Polycarbonate with molded-in graduations. For low-speed centrifugation up to 6000 x g. Autoclavable. The 15 mL size is used with Settlometer Kit (P/N 21157-00). Supplied in packages of 10, without closures. Nalge.

Technické údaje

Délka:	120 mm
Kapacita:	15 mL
Obsah balení:	10 /pk