



Test kit, ascorbic acid, model ASC-1

Kat. číslo produktu:

2308100

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno



Model ASC-1 Ascorbic Acid Test Kit was developed in cooperation with beverage producers and uses the drop-count titration method to measure ascorbic acid (vitamin C). Utilizes drop count titration with a range of 10-200 mg/L. Contains approximately 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Metoda:	Drop count titration/ Iodometric
Parametr:	Ascorbic Acid
Platforma:	Drop Count Titration
Počet testů:	100
Rozsah měření:	10 - 200 mg/L
Smallest increment:	10
Styl kufříku:	D
Typ testovací soupravy:	ASC-1
Vysvětlivka:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.