



Test kit, molybdate, model MO-LR

Kat. číslo produktu:

2359300

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno



Molybdates are gaining over chromates for use as corrosion inhibitors due to lower toxicity. The high cost of these treatments warrant careful control of their concentration. Kit reagents react with Molybdate to form a yellow color, which is measured with a color disc.

Kit reagents react with Molybdate to form a yellow color, which is measured with an easy-to-use color disc method

Effectively monitors corrosion inhibitors

Technické údaje

Metoda:	Color disc/ Ternary complex
Parameter:	Molybdate
Parametr:	Molybdate - As Molybdenum
Platforma:	Barevný disk
Počet testů:	100
Přístroj:	MO-LR
Rozsah měření:	Single Range
Smallest increment:	0.25
Standardní rozsah měření:	0 - 3.0 mg/L Mo
Styl kufříku:	D
Typ testovací soupravy:	MO-LR