



Agar LES, m-Endo, dehydrated 500 g

Kat. číslo produktu:

2406634

Nebezpečné

ZASTARALÁ POLOŽKA

Tato položka již není dostupná.



Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

Technické údaje

Media Type:	m-Endo LES Agar
Obsah balení:	500 g