



Filter, TSS glass fiber, pore size 1.5 µm, diameter 55 mm, 100/pk

Kat. číslo produktu:	253055
CZK Cena bez DPH:	Kontaktujte nás
Dodání do 2 týdnů	

TSS Filters are recommended for use in Standard Methods 2540 C and 2540 D. Optimized surface characteristics of this binderless borosilicate glass microfiber filter eliminates pin-holing and fiber-shedding to ensure consistent, accurate results. Fast flow rate, fine porosity and 1.5 µm particle retention makes the TSS Filter the preferred choice of environmental chemists for total suspended solids analysis of water and wastewater.

Technické údaje

Materiál:	Glass Microfiber
Obsah balení:	100 /pk
Průměr:	55 mm
Velikost póru:	1.5 µm