



Sample Container, 120 mL High PF

Kat. číslo produktu:	2552200
CZK Cena bez DPH:	Kontaktujte nás
Dodání do 2 týdnů	

These 120 mL (4oz.) High Profile Corning Snap-Seal™ polypropylene containers are designed for a wide variety of applications and provide a reliable leak-tight seal when closed properly. Approximately 3.5 inches tall (91 mm) and 1.75 inches (45 mm) diameter. Plastic, snap seal, translucent. Sold individually.

Technické údaje

Kapacita:	120 mL - 1 each
Obsah balení:	1 /pk
Průměr:	45 mm (approx)