



Sample container, 120 mL, pk/250

Kat. číslo produktu:	2552212
CZK Cena bez DPH:	Kontaktujte nás
Dodání do 2 týdnů	

Plastic, leak-proof. Corning Snap Seal offers a non-sterile, leak-proof, attached cap container for field use in sampling water, soils, and sludge. Locking device helps prevent accidental opening of sample; container closes with patented sealing mechanism so sample won't leak.

Technické údaje

Kapacita:	120 mL - 250/pk
Obsah balení:	250 /pk
Průměr:	45 mm (approx)