



Cylinder, graduated, 1000 mL, w/ handles

Kat. číslo produktu:	2612953
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

One-piece polypropylene construction with a molded in handle top ensure a secure, solid grip and superior control for safer, more convenient pouring. The “Total Contained, Total Delivered ”double scale is accurate and easy to read (cylinder will not create a meniscus), while the hexagonal base prevents tipping or rolling. Hand opening is 89x 29 mm (3 x 11/8"). Steam autoclavable upto 121°C (250°F).

Technické údaje

Certifikace:	None
Kapacita:	1000 mL
Materiál skříně:	Polypropylene
Metoda měření:	To Deliver/To Contain
Obsah balení:	1 each
Typ:	S dvojitou stupnicí