



Bottle, wash, diamond alert, 500 mL, Variety 6/pk

Kat. číslo produktu:	2613049
CZK Cena bez DPH:	Kontaktujte nás
Dodání do 2 týdnů	

Wide-mouthed and color-coded wash bottles provide safe and positive identification of common laboratory solutions. Standard Diamond Alert symbols indicate type of hazard contents may represent. Polyethylene bottles with color-coded polypropylene screw caps and spray-tip assembly. Variety 6-pack includes one bottle each for water, acetone, methanol, isopropanol, toluene, and ethanol.

Technické údaje

Barva:	Transparent
Jednotka:	6 pk
Kapacita:	500 mL
Materiál:	LDPE
Obsah balení:	6 pk
Stem:	Blue Cap/Stem