



## Filter, membrane, pore size 0.2 µm, diameter 47 mm, 200/pk (Gelman)

**Kat. číslo produktu:** 2618800  
**CZK Cena bez DPH:** Kontaktujte nás  
**Datum odeslání neuvedeno**

For isolation and enumeration of stressed organisms (*Pseudomonas* species). High retention and flow rate. Certified for use in membrane filtration analysis. White polyethersulfone media. STerile and individually packaged. Pore size: 0.2 µm. Diameter: 47mm.

---

### Technické údaje

|                |                  |
|----------------|------------------|
| Filter Type:   | Filter           |
| Materiál:      | Polyethersulfone |
| Obsah balení:  | 200 /pk          |
| Průměr:        | 47 mm            |
| Velikost póru: | 0.2 µm           |