



Thermometer, Enviro-Safe, Non-Mercury, -10 - 225 °C

Kat. číslo produktu:	2635700
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

Enviro-Safe non-Mercury thermometers feature an exclusive biodegradable liquid and dye. Easy to read, aesthetically distinctive due to dark green liquid and permanently fused markings. Each thermometer, individually calibrated against NIST traceable equipment, identified by unique serial number for ISO

9000 or other quality systems.

Technické údaje

Dělení teploměru:	1 °C
Délka:	405 mm
Obsah balení:	1 each
Teplotní rozsah:	-10 °C - 225 °C
Typ kapalinového teploměru:	Non-Mercury