



Cylinder, Graduated, 1000 mL, Certified

Kat. číslo produktu:	2636253
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

Certified Class A, this glass cylinder is calibrated "to deliver" according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in white marking. Bumper guard included for all sizes ≥ 25 mL.

Technické údaje

Certifikace:	Class A
Kapacita:	1000 mL
Materiál skříně:	Glass
Metoda měření:	To Deliver
Obsah balení:	1 each
Typ:	S jednou stupnicí
Výška:	465 mm