



Cylinder, graduated, mixing, certified, 25 mL

Kat. číslo produktu:	2636340
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

This certified, Class A glass mixing cylinder is calibrated “to contain” according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in durable white markings with funnel-top. Glass penny head stopper included.

Technické údaje

Certifikace:	Class A
Kapacita:	25 mL
Materiál skříně:	Glass
Metoda měření:	To Contain
Obsah balení:	1 each
Typ:	S jednou stupnicí, směšovací