



Iron (Ferrous) Color Disc Test Kit, Model IR-18C, 0-7 mg/L

Kat. číslo produktu:

2667200

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Dodání do 2 týdnů



The most accurate visual test kit method, Hach color discs provide fast, accurate comparisons. Range: 0-7 mg/L. Contains 100 tests.

Simple, fast, color disc method

Rugged test kit contains all required reagents and apparatus to perform analysis in the field

Technické údaje

Metoda:	Color disc/ 1,10 Phenanthroline
Parameter:	Iron Ferrous
Parametr:	Iron, medium range - As Fe ²⁺ (ferrous iron)
Platforma:	Barevný disk
Počet testů:	100
Přístroj:	IR-18C
Rozsah měření:	High Single Range
Smallest increment:	0.2
Styl kufříku:	B
Typ:	S jednou stupnicí
Typ testovací soupravy:	IR-18C