



Test kit, RapidSilver[™] visual

Kat. číslo produktu: 2674500

ZASTARALÁ POLOŽKA
Tato položka již není dostupná.

Test Kit Silver, Model RapidSilver[™]. Range 0-50 ppb. Silver As Ag Visual Comparison.

Provides fast, accurate information on silver discharges

Used by surface finishers, photographic film manufacturers and processors, and others to prevent pollution and economic loss

Kit directions also include digestion and dilution procedures needed for testing silver strike bath samples

Technické údaje

Metoda:	Visual Comparison
Parameter:	Silver
Parametr:	Silver-As Ag
Počet testů:	100
Přístroj:	NA
Rozsah měření:	5 - 50 ppb Ag
Smallest increment:	5 µg/L
Standardní rozsah měření:	5 - 50 ppb Ag
Styl kufříku:	F
Typ testovací soupravy:	RapidSilver [™]