



Test kit, hypochlorite, model CN-HRDT

Kat. číslo produktu:

2687100

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno



Hypochlorite Kit, Model CN-HRDT with Digital Titrator

Utilizing the Hach Digital Titrator, this kit is used to determine hypochlorite, high range (Bleach) solution concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 5 - 15% as Cl₂: 50 - 150 g/L as Cl₂. Contains 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Designed for on-site analysis of bulk hypochlorite deliveries or bulk storage solutions at drinking water or chloramination facilities

Verify concentration of hypochlorite solutions at time of delivery

Kit contains all required reagents, solutions, and glassware

Technické údaje

Číslo metody:	10100
Metoda:	Digital Titrator/ Thiosulfate
Parametr:	Hypochlorite, high range (Bleach) - As Cl ₂
Platforma:	Digital Titrator
Počet testů:	100
Rozsah měření:	5 - 15 % = 50-150 g/L Cl ₂
Smallest increment:	0.05%
Styl kufříku:	G
Typ testovací soupravy:	CN-HRDT