



## Plastic ampules, pseudomonas broth, pk/50

|                      |                 |
|----------------------|-----------------|
| Kat. číslo produktu: | 2812250         |
| CZK Cena bez DPH:    | Kontaktujte nás |
| Dodání do 2 týdnů    |                 |

Pseudomonas broth plastic ampules, 2 mL. For detection of Pseudomonas species. pk/50.

### Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

### Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

### Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

---

## Technické údaje

|               |                   |
|---------------|-------------------|
| Materiál:     | Plastic           |
| Media Type:   | Pseudomonas Broth |
| Obsah balení: | 50 50/pk          |