



Plastic ampules R2A, pk/50

Kat. číslo produktu: 2812350
CZK Cena bez DPH: Kontaktujte nás
Datum odeslání neuvedeno

R2A Broth Ampules; Plastic

R2A Plastic Ampules for detection of stressed heterotrophic bacteria; pk/50.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of

your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make

wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking

water samples, Hach has prepared media and convenient apparatus to make

testing easy and fast.

Technické údaje

Materiál:	Plastic
Media Type:	R2A
Obsah balení:	50 /pk