



Vacuum filter microfunnel, pk/50.

Kat. číslo produktu:	2831500
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

Feature a completely assembled filtration unit: Funnel, filter, nutrient pad and petri dish. Filter: GN-6 Metrical with 0.45 µm pore size and 47 mm diameter; white, gridded surface. Funnel: Polypropylene Petri Dish: polystyrene lid and polypropylene base Nutrient Pad: Cellulose. Fits 2831500, 2846100 and 2846150 manifolds. Sterile and individually bagged.

Technické údaje

Materiál:	Mixed-Cellulose Ester
Obsah balení:	50
Průměr:	47 mm
Velikost póru:	0.45 µm