



Bottle, BOD, 300 mL, 24/pk, serialized: 1-24

Kat. číslo produktu:	2898700
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included Serialized bottles have silk-screened, sequential numbers (1-24) and include a bar code.

Technické údaje

Body Material:	Glass
Kapacita:	300 mL
Materiál:	Borosilicate Glass
Obsah balení:	24 /pk
Package size UOM:	24/pk, Serialized: 1 - 24