



pH paper, 1.0 - 11.0 pH range

Kat. číslo produktu:	39133
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

Provides simple, fast, economical method of pH measurement. Sharp color change for each 0.5 pH unit. Convenient for sample pH adjustments. Recommended for sewage water, not finished water. Five 4.6-meter (15') rolls per package.

Technické údaje

Obsah balení:	5 4.6-meter (15') rolls per package
Rozsah měření:	1.0 - 11.0 pH units