



Beaker, Griffin, Low Form, Glass, 600 mL

Kat. číslo produktu:	50052
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

Technické údaje

Jednotky:	1 each
Kapacita:	600 mL
Materiál:	Glass, Griffin, Low Form
Materiál skříně:	Glass