



## QL filters, pore size 8 - 12 $\mu\text{m}$ , diameter 70 mm, circular, 100/ pk, qualitative

Kat. číslo produktu: 50653

**ZASTARALÁ POLOŽKA**  
Tato položka již není dostupná.

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

---

### Technické údaje

|               |           |
|---------------|-----------|
| Materiál:     | Cellulose |
| Obsah balení: | 100 /pk   |
| Průměr:       | 70 mm     |