



## **Xenomax TSS-FVS Filter 47 mm, 100/pk**

<b>Kat. číslo produktu:</b>	<b>2511307</b>
CZK Cena bez DPH:	Kontaktujte nás
Dodání do 2 týdnů	

Xenomax<sup>™</sup> TSS-FVS Filters are specifically designed to ensure compliance with the filter preparation requirements of Standard Methods 2540 D and Standard Methods 2540 E (Fixed and Volatile Solids Ignited at 550°C). Independent 3rd party evaluation of Xenomax and 934-AH<sup>®</sup> filters in both methods give demonstrably equivalent results.

In addition, these 1.5 µm binder-free, fast flowing, fine porosity borosilicate glass microfiber filters have been optimized for better wet strength/handleability and high temperature (550°C) blank stability. FVS Filters are available in popular diameter sizes and ready for use in a wide variety of labware including Buchner funnels and Gooch Crucibles.

934-AH is a registered trademark of Whatman, Inc./GE Healthcare.

---

### **Technické údaje**

Materiál:	Glass Microfiber
Obsah balení:	100 /pk
Průměr:	47 mm
Velikost póru:	1.5 µm