



Be Right™



De-bubbler including tubing and pressure relief valves for ScrubberMonitor cabinets

Kat. číslo produktu: **APPS10011888**

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno

Complete de-bubbler kit for the ScrubberMonitor EGCS washwater monitoring cabinets. Made with marine-grade PVC built to withstand harsh environments, this system is designed to take air out of the incoming sample to provide consistent sample flow to the ScrubberMonitor. Includes all of the piping, pressure regulators, and flexible tubing needed to connect to the cabinet.