



## Reagent Kit for EZ4041 Total Hardness Analyser (acid, buffer, colour and EDTA solution), 10-250 mg/L, 10-25% range

**Kat. číslo produktu:**

**APPC4041KTOL-EN** ~~Nebezpečné~~

**CZK Cena bez DPH:**

Kontaktujte nás

**Datum odeslání neuvedeno**



### Ready-to-use Reagents for EZ Series Online Analysers

#### Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

#### Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

#### Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

---

### Technické údaje

Jaký je obsah balení:

APPC4041-01-EU Acid Solution

APPC4041-02-EU Buffer Solution

APPC4041-03 Colour Solution

APPC4041-04 EDTA Solution (0.005 M)

Parametr:

Total Hardness

---

### Jaký je obsah balení

APPC4041-01-EU Acid SolutionAPPC4041-02-EU Buffer SolutionAPPC4041-03 Colour SolutionAPPC4041-04 EDTA Solution (0.005 M)