




Reagent Kit for EZ4043 Total Hardness Analyser (acid, buffer, colour and EDTA solution), 0.25-10 mg/L CaCO₃

Kat. číslo produktu:

APPC4043KTO-EU **Nebezpečné**

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno



Ready-to-use Reagents for EZ Series Online Analysers

Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

Technické údaje

Jaký je obsah balení:

APPC4041-01-EU Acid Solution

APPC4041-02-EU Buffer Solution

APPC4041-03 Colour Solution

APPC4043-04 EDTA Solution, 0.005 N

Parametr:

Total Hardness

Jaký je obsah balení

APPC4041-01-EU Acid SolutionAPPC4041-02-EU Buffer SolutionAPPC4041-03 Colour SolutionAPPC4043-04 EDTA Solution, 0.005 N